



# Superheating Mizorogi-Heck reaction: 2-liquid mixture reaction



## Overview

The basic 2-line system performs library synthesis under superheating in this application example. The **OptimFlow** allows safe and easy determination of reaction conditions, such as the concentrations of reaction solutions, equivalents of catalysts, and reaction temperatures.



## Preparation of Reagents

**Reagent A:** Iodobenzene **1a** (50.0 mg, 0.245 mmol) was dissolved in 1 mL of MeCN to have a concentration of 0.24 M. [**1b**: 1-chloro-4-iodobenzene (56.1 mg, 0.245 mmol), **1c**: 4-iodobenzonitrile (58.4 mg, 0.245 mmol), **1d**: 1-iodo-4-methylbenzene (53.4 mg, 0.245 mmol), **1e**: 1-iodo-4-methoxybenzene (57.4 mg, 0.245 mmol), and **1f**: 1-iodo-4-(trifluoromethyl)benzene (66.7 mg, 0.245 mmol) ]

**Reagent B:** Butyl acrylate **2** (47.1 mg, 0.368 mmol), N,N-disopropylethylamine (47.5 mg, 0.368 mmol), and palladium(II) acetate (0.275 mg, 1.23  $\mu$ mol) were dissolved in 1 mL of MeCN.

## Device Setup

Line 1, System C ver. 3.1; Line 2, System C ver. 3.1

(Our compact back pressure valve was used for the plunger pump Out to improve the constant flow rate.  
Discharge pressure: 1.7 MPa)

From Loop to Reactor: A SUS tube with a length of 500 mm and an I.D. of 0.50 mm

BPR = 0.7 MPa, the 2-reagent mode with a T-shaped mixer

The BPR was set to 0.7 MPa. Liquid feeding started at room temperature, which was then raised to 160 °C in an oil bath.

The experiment started when the flow velocity stabilized.

Batch synthesis has  
an upper limit of 82 °C  
for heating with reflux.

## Reactors

A SUS tube (I.D. 1.0 mm, volume 2 mL) was placed in an oil bath and connected to the **OptimFlow** as a reactor. A SUS tube (500 mm) was connected immediately after the reactor for air cooling. A PTFE tube (500 mm) was used for the connection to ABPR (From Reactor 1 to BPR: a length of 1000 mm and an I.D. of 0.50 mm).

## Fraction Collector

In each experiment, the reaction solutions were collected in a test tube, with the Pre Stream Ratio set to 800  $\mu$ L, the Post Stream Ratio to 800  $\mu$ L, and the Fraction Volume to 2000  $\mu$ L.

Experiment

The parameters for each experiment were set as shown in the software input example (Experimental parameter). The reference amount of Reagent A (Volume of Reagent 1) was 200  $\mu\text{L}$ , and the Residence Time was 2.5 min. The flow rates, the amounts of reagents used, and the total amount of reaction solutions in the above settings can be viewed in the Calculated value tab.

Experimental parameter	Calculated value										Procedure & details	Common parameter
	Expt.1	Expt.2	Expt.3	Expt.4	Expt.5	Expt.6	Expt.7	Expt.8	Expt.9	Expt.10		
Reagent1	1A1	1A2	1A3	1A4	1A5	1A6	1A7	1A8	1A9	1A10		
Reagent2	2A1	2A1	2A1	2A1	2A1	2A1	2A1	2A1	2A1	2A1		
Residence Time (min)	2.5000	2.5000	2.5000	2.5000	2.5000	2.5000	2.5000	2.5000	2.5000	2.5000		
Volume Of Reagent1 ( $\mu\text{L}$ )	200	200	200	200	200	200	200	200	200	200		
Conc.of Reagent1 (M)	0.24	0.24	0.24	0.24	0.24	0.24	0.24	0.24	0.24	0.24		
Conc.of Reagent2 (M)	0.24	0.24	0.24	0.24	0.24	0.24	0.24	0.24	0.24	0.24		
Mol Ratio of Reagent1	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00		
Mol Ratio of Reagent2	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00		
Pre Stream ( $\mu\text{L}$ )	800	800	800	800	800	800	800	800	800	800		
Post Stream ( $\mu\text{L}$ )	800	800	800	800	800	800	800	800	800	800		
Fraction Volume ( $\mu\text{L}$ )	2000	2000	2000	2000	2000	2000	2000	2000	2000	2000		
Set BPR (MPa)	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7		
Upper Limit of Pressure (MPa)	3.00	3.00	3.00	3.00	3.00	3.00	3.00	3.00	3.00	3.00		

After the reactions, 5  $\mu\text{L}$  of each resulting reaction solution was diluted with 120  $\mu\text{L}$  of MeOH and analyzed with a UPLC-MS. The yields were calculated using area normalization with UV areas of the products, raw materials, and byproducts.

Entry	1	R	3 (%) <sup>*</sup>	1 (%) <sup>*</sup>
1	1a	H	86	14
2	1b	Cl	79	21
3	1c	CN	75	11
4	1d	Me	74	26
5	1e	OMe	55	45
6	1f	CF <sub>3</sub>	61	34

\*The yields were analyzed and calculated for a sample of 1:1 molar ratio mixture of 3 and 1 at wavelengths with equal UV intensity ratios of 3 and 1.

Analysis wavelengths (nm)


R=H: 220, Cl: 225, CN: 262, Me: 225, OMe: 228, CF<sub>3</sub>: 244

Analysis conditions


Column: ACQUITY UPLC BEH C18(1.7  $\mu\text{m}$ , 2.1  $\times$  30mm), 60  $^{\circ}\text{C}$   
Solvent A: 0.05% TFA aq., B: 0.05% TFA/MeCN  
Gradient: 30-100% B in 2.2 min Flow Rate: 1.0 mL/min

An example of condition determination for R = Cl is shown below. <Line 1: 1b, 2, DIEA, Line 2: Pd(OAc)<sub>2</sub>>


(i) Relationship between reaction concentration and catalyst concentration



(ii) Relationship between catalyst concentration and residence time



(iii) Reaction temperature



✓ The **OptimFlow** is designed for multiple samples with liquid handlers and thus requires no reparation of reaction solutions for determining the reaction concentration and the catalyst concentration.

✓ The residence time can be changed only by input on the PC.